

VECTOR BORNE DISEASE

STANDARD Q Chikungunya IgM/IgG



Test type	Professional Use Only
Intended use	Detection of Chikungunya IgM and IgG antibodies
Specimen type	Whole blood, Serum, Plasma
Specimen volume	10µl
Testing time	15-20 mins
Storage condition	2 – 40°C / 36 - 104°F

» Test Performance

Reference : Internal evaluation

	Reference	Sensitivity	Specificity
IgM	ELISA	100% (21/21)	97.7% (253/259)
IgG	ELISA	100% (21/21)	99.6% (258/259)

» Ordering Information

Product	Pack Size	CAT No.
Q Chikungunya IgM/IgG Test	25 Tests	09CHI20D