

STANDARD i-Q COVID-19 Ag Home Test

Test type	Self-diagnostic test
Specimen type	Nasal swab
Testing time	15-30 mins
Storage condition	2-30°C / 36-86°F



» Ordering Information

Product	Pack Size	CAT No.
i-Q COVID-19 Ag Home Test	1 test	09COV121
i-Q COVID-19 Ag Home Test	2 test	09COV120H
COVID-19 Ag Control Swab	Pos x 10, Neg x 10	10COVC11

STANDARD Q Influenza A/B



Test type	Professional Use Only
Intended use	Detection of influenza A/B antigens
Specimen type	Nasopharyngeal swab / Nasopharyngeal wash/aspirate
Testing time	8-12 mins (Do not read after 30 mins)
Storage condition	2 – 40°C / 36 - 104°F



» Test Performance

Reference : Internal evaluation

Reference	Sensitivity	Specificity
PCR	A: 97.44% (95% CI: 86.52-99.94%)	A: 100% (95% CI: 99.12-100.00%)
	B: 90.63% (95% CI: 74.98-98.02%)	B: 98.82% (95% CI: 97.26-99.61%)

» Ordering Information

Product	Pack Size	CAT No.
Q Influenza A/B Test	25 Tests	09INF40D
Influenza A/B Control	Pos x 10 / Neg x 10	10INFC10

STANDARD Q RSV Ag



Test type	Professional Use Only
Intended use	Detection of Respiratory Syncytial Virus (RSV) antigens
Specimen type	Nasopharyngeal swab / Nasopharyngeal wash/aspirate
Testing time	15 mins (Do not read after 30 mins)
Storage condition	2 – 40°C / 36 - 104°F



» Test Performance

Reference : Internal evaluation

Reference	Sensitivity	Specificity
PCR	92.45% (49/53)	98.44% (126/128)

» Ordering Information

Product	Pack Size	CAT No.
Q RSV Ag	25 Tests	09RSV40D
RSV Ag Control	Pos x 10 / Neg x 10	10RSVC10